



Support Your  
Healthy Lifestyle



# Compression Stockings

Relax the Leg and Help Improve Blood Circulation



**USA**  
**OPPO MEDICAL INC.**  
Seattle City Center  
1420 Fifth Ave., Ste. 220080,  
Seattle, WA 98101, USA  
Tel: 1-206-575-8843  
Fax: 1-206-575-9351  
Website: [www.oppomedical.com](http://www.oppomedical.com)  
E-mail: [sales@oppomedicalusa.com](mailto:sales@oppomedicalusa.com)

**SOUTH AMERICA**  
**OPPO MEDICAL S.R.L.**  
Avenida Artigas 2164,  
Asunción, Paraguay  
Tel: 595-21-282368  
Fax: 595-21-291430  
E-mail: [anva@pla.net.py](mailto:anva@pla.net.py)

**CHINA**  
**PLUSON (SHANGHAI) CO. LTD.**  
15A, Huamin Empire Plaza,  
726 Yan An West Rd.,  
Shanghai 200050, China  
Tel: 86-21-52391585  
Fax: 86-21-52395055  
E-mail: [sales@oppomedical.com.cn](mailto:sales@oppomedical.com.cn)

**TAIWAN**  
**OPPO MEDICAL CORP.**  
9F, No. 297, Sec.4, Chung Hsiao E. Rd.,  
Taipei 10696, Taiwan  
Tel: 886-2-27788710  
Fax: 886-2-27739935  
E-mail: [sales@oppomedical.com](mailto:sales@oppomedical.com)

OPPO and OPPO logo are trademarks of  
OPPO MEDICAL INC. USA © 2011





## About OPPO

Established in Seattle, Washington in 1997, OPPO Medical has built a powerful network of affiliate companies in America and Asia, and has more than 70 distributors worldwide. We offer the most extensive product line of functional supports with:

- + Long standing commitment to uncompromising quality in workmanship and customer service.
- + Keep understanding of the rehabilitation market and the needs of our customers.
- + Passion for innovation, shared by a team of highly-motivated and talented employees.
- + Powerful network of affiliate companies and factories in the US, Taiwan, China and Paraguay, positioned to serve you, our distributors, whether you are in North America, South America, Central America, Europe, Asia, Middle East or Africa.

We continue to use technology to streamline manufacturing processes and to foster networking communication for our distributors worldwide. We are expanding the frontiers of orthopedic medicine with innovative materials and designs, and a holistic approach to wellness.

## About OPPO Compression Stockings

OPPO Medical is committed to bringing you high-quality, attractive, reliable vein protection at a value price. OPPO Compression Stockings help improve blood circulation and reduce symptoms associated with vein disorders. Made for both men and women; we offer a four series line of stockings - **Caring Series**, **Voyage Series**, **Eminence Series** and **Anti-embolism Series** - for ideal care and management of venous diseases to meet the needs of a wide range of customers. With a contoured design, soft texture and various compression levels, OPPO Compression Stockings ensure a perfect fit and optimal compression for greater comfort.



## CONTENT

OPPO Compression Stockings Features .....	3
Graduated Compression Benefits .....	5
Caring Series .....	8
Voyage Series .....	14
Eminence Series .....	16
Anti-embolism Series .....	19
Appendix .....	21







Support Your  
Healthy Lifestyle

## OPPO Compression Stockings Features

OPPO Compression Stockings utilize extraordinary LYCRA® fiber and German-engineered professional knitting machines to provide ideal management of vein problems and prevent symptoms. LYCRA® is high-quality, highly-resilient elastic material that recovers well after use for superior performance, durability and comfort. OPPO Compression Stockings provide stable gradient pressure and precise pressure stages for different parts on lower limbs.

## Various Designs for Maximum Comfort

OPPO Compression Stockings' comfortable texture and various thoughtful designs are suitable for daily use for long periods of time and may be work during most daily activities, e.g. working, traveling, or special medical requirements of vein related conditions, consistent with physician's recommendation.



### Pantyhose

OPPO Pantyhose offers an ergonomically designed, low-front, high-back waistband that comfortably and perfectly fits the hips.



### Pantyhose

The high-ventilation mesh gusset helps ensure all-day comfort.



### Stockings with lace top band

Lace top with silicone strips prevent slipping and binding.



### Stockings with silicone top band

The stocking and top band are sewn together by fine, flat seams to maximize comfort and reduce chafing. Wide band with silicone strip provides a soft touch to the skin while maintaining secure fit to prevent slipping.




### All Stockings

The reciprocated knit on heel and toes ensures extra durability. The contoured design of heel space provides better fit and comfort without wrinkling.



## Graduated Compression Benefits

If you are standing or sitting for long periods of time, the calf muscles are not contracting to squeeze the veins carrying blood from the lower limbs back to the heart. OPPO Compression Stockings help promote blood circulation and improve valve function by applying graduated compression on the lower limbs. The pressure is strongest at the ankle and gradually decreases up the leg. The gradient compression exerted from the stockings provides a pumping effect to accelerate blood flow in the legs upward to the heart, reducing the accumulation of blood in the lower limbs and eliminating resultant edema.



## Compression Level and Indications

OPPO Compression Stockings provide different strengths of compression levels from Class 1 to Class 3. The classes are prescribed in millimeters of mercury (mmHg), which represents the amount of pressure put on the legs by the stockings. Proper fit and pressure are essential for compression stockings to work effectively, consumers should follow doctor's recommendation for the correct level of compression.



- Edema / Swollen legs
- Muscle weakness
- Leg pain
- Mild varicosity
- Prevent varicosity during pregnancy
- Prevent Economy Class Syndrome



- Edema / Swollen legs
- Spider veins
- Moderate varicose veins
- Post-surgical wear
- Varicose veins during pregnancy
- Economy Class Syndrome



- Severe edema
- Scleroderma
- Recurrent varicose veins
- Chronic venous insufficiency
- Post-surgical wear
- Severe varicosity during pregnancy

## Varicosity and Its Symptoms

Varicosity results from a chronic increase in venous blood pressure, which dilates the vein. When the valves in the vein cannot function properly, blood flow reverses. The blood accumulates in the vein, increasing venous pressure and the possibility of congestion while causing the vein to bulge and twist, i.e. varicosity.

The most frequent signs of varicosity:

### 1. Cosmetic symptoms:

- + Discolored skin
- + Peeling skin
- + Skin ulcers
- + Superficial vein become visible

### 2. Physical symptoms:

- + Ankle edema
- + Muscle weakness
- + Itching or aching of the leg

### 3. Functional symptoms:


- + Heaviness of the leg
- + Cannot sit or stand for a long time
- + Irritation in prolonged positions

## Compression Therapy –

The most fundamental treatment and preventive measure for venous disorders


### Healthy Vein

Normal valve allows blood to flow towards the top and the closed venous valves prevent the reverse flow of the blood.




### Dysfunctional Vein

The venous valves cannot close effectively. The blood flows back and clots in the vein, which distends the section of vein below it, causing further valves to fail.



### Dysfunctional Vein with OPPO Compression Stockings

Compression stockings compress the veins improving the flow of blood.





# Reading Instructions



2830

**KNEE HIGH COMPRESSION STOCKINGS**

## Indication

- Mild varicosity
- Leg pain, edema and muscle weakness
- Pregnancy

Size: I - VII

Color 



## Product Number & Name

Provide quick search for specific product.

## Indications

## Size Range

## Color


## Compression Level

Compression level of each product. (from Class 1 to Class 3)

## Toe Design

Open-toe or Closed-toe.

Product Image



## SIZING:

Measure circumference of the ankle, right above the ankle bones (malleoli) (cB) and around the fullest part of the calf (cC). Measure knee height (ID).

SIZE	I	II	III	IV	V	VI	VII
cC	31 - 35 cm 12 1/4 - 13 3/4"	33 - 38 cm 13 - 15"	35 - 40 cm 13 3/4 - 15 3/4"	37 - 42 cm 14 1/2 - 16 1/2"	39 - 44 cm 15 1/4 - 17 1/4"	41 - 47 cm 16 1/8 - 18 1/2"	43 - 49 cm 16 7/8 - 19 1/4"
cB	18 - 20 cm 7 - 7 7/8"	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"
ID	35 - 39 cm 13 3/4 - 15 3/8"	35 - 39 cm 13 3/4 - 15 3/8"	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	40 - 45 cm 15 3/4 - 17 5/8"	40 - 45 cm 15 3/4 - 17 5/8"

## Measurements Chart

Precise measuring of circumference and length is essential to having a correct fitting of stockings. The measurements chart shows the few measurements that are needed to determine the right size for each type of compression stockings.

## Caring Series

OPPO Caring Compression Stockings provide a wide selection of designs with different compression levels and lengths for individual fit for daily activities or during pregnancy. The ergonomically contoured knit is specially designed for women to provide precise fit and optimum compression for the legs to accelerate blood circulation, reduce varicose vein, relieve foot fatigue and help prevent swelling conditions in all ages.





2830

**KNEE HIGH COMPRESSION STOCKINGS**

**Indication**

- Mild varicosity
- Leg pain, edema and muscle weakness
- Pregnancy

Size: I - VII

Color



2824

**KNEE HIGH COMPRESSION STOCKINGS**

**Indication**

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear

Size: I - VII

Color



2802

**KNEE HIGH COMPRESSION STOCKINGS**

**Indication**

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear

Size: I - VII

Color



2803

**KNEE HIGH COMPRESSION STOCKINGS**

**Indication**

- Scleroderma
- Severe edema
- Recurrent varicose veins
- Severe varicose veins during pregnancy
- Post-surgical wear
- Chronic venous insufficiency

Size: I - VII

Color



2852

**MID THIGH COMPRESSION STOCKINGS**

**Indication**

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear

Size: I - VII

Color



2851

**MID THIGH COMPRESSION STOCKINGS**

**Indication**

- Mild varicosity
- Leg pain, edema and muscle weakness
- Pregnancy

Size: I - VII

Color



2840


**MID THIGH COMPRESSION STOCKINGS**

**Indication**

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear

Size: I - VII


Color



**SIZING:**

Measure circumference of the ankle, right above the ankle bones (malleoli) (cB) and around the fullest part of the calf (cC). Measure knee height (ID).

SIZE	I	II	III	IV	V	VI	VII
cC	31 - 35 cm 12 1/4 - 13 3/4"	33 - 38 cm 13 - 15"	35 - 40 cm 13 3/4 - 15 3/4"	37 - 42 cm 14 1/2 - 16 1/2"	39 - 44 cm 15 1/4 - 17 1/4"	41 - 47 cm 16 1/8 - 18 1/2"	43 - 49 cm 16 7/8 - 19 1/4"
cB	18 - 20 cm 7 - 7 7/8"	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"
ID	35 - 39 cm 13 3/4 - 15 3/8"	35 - 39 cm 13 3/4 - 15 3/8"	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	40 - 45 cm 15 3/4 - 17 5/8"	40 - 45 cm 15 3/4 - 17 5/8"



**SIZING:**

Measure circumference of the ankle, right above the ankle bones (malleoli) (cB) and around the fullest part of the calf (cC). Measure knee height (ID).

SIZE	I	II	III	IV	V	VI	VII
cC	31 - 35 cm 12 1/4 - 13 3/4"	33 - 38 cm 13 - 15"	35 - 40 cm 13 3/4 - 15 3/4"	37 - 42 cm 14 1/2 - 16 1/2"	39 - 44 cm 15 1/4 - 17 1/4"	41 - 47 cm 16 1/8 - 18 1/2"	43 - 49 cm 16 7/8 - 19 1/4"
cB	18 - 20 cm 7 - 7 7/8"	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"
ID	35 - 39 cm 13 3/4 - 15 3/8"	35 - 39 cm 13 3/4 - 15 3/8"	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	40 - 45 cm 15 3/4 - 17 5/8"	40 - 45 cm 15 3/4 - 17 5/8"





2871

**THIGH HIGH COMPRESSION STOCKINGS**

**Indication**

- Mild varicosity
- Leg pain, edema and muscle weakness
- Pregnancy

Size: I - VII

Color 



2865

**THIGH HIGH COMPRESSION STOCKINGS**

**Indication**

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear

Size: I - VII

Color 



2891

**COMPRESSION PANTYHOSE**

**Indication**

- Mild varicosity
- Leg pain, edema and muscle weakness

Size: I - VII

Color 



2880

**COMPRESSION PANTYHOSE**

**Indication**

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear

Size: I - VII

Color 



2875

**THIGH HIGH COMPRESSION STOCKINGS**

**Indication**

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear

Size: I - VII

Color 




2861

**THIGH HIGH COMPRESSION STOCKINGS**

**Indication**

- Scleroderma
- Severe edema
- Recurrent varicose veins
- Severe varicose veins during pregnancy
- Post-surgical wear
- Chronic venous insufficiency

Size: I - VII

Color 



2882

**COMPRESSION PANTYHOSE**

**Indication**

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear

Size: I - VII

Color 



2884


**COMPRESSION PANTYHOSE**

**Indication**

- Scleroderma
- Severe edema
- Recurrent varicose
- Post-surgical wear
- Chronic venous insufficiency

Size: I - VII


Color 



**SIZING:**

Measure circumference of the ankle, right above the ankle bones (malleoli) (cB), around the fullest part of the calf (cC) and around the fullest part of the thigh (cG). Measure thigh height (IG).

SIZE	I	II	III	IV	V	VI	VII
cG	45 - 53 cm 17 5/8 - 20 7/8"	48 - 56 cm 18 7/8 - 22"	51 - 59 cm 20 - 23 1/8"	54 - 63 cm 21 1/4 - 24 3/4"	57 - 66 cm 22 3/8 - 26"	60 - 70 cm 23 5/8 - 27 1/2"	63 - 73 cm 24 3/4 - 28 3/4"
cC	31 - 35 cm 12 1/4 - 13 3/4"	33 - 38 cm 13 - 15"	35 - 40 cm 13 3/4 - 15 3/4"	37 - 42 cm 14 1/2 - 16 1/2"	39 - 44 cm 15 1/4 - 17 1/4"	41 - 47 cm 16 1/8 - 18 1/2"	43 - 49 cm 16 7/8 - 19 1/4"
cB	18 - 20 cm 7 - 7 7/8"	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"
IG	62 - 70 cm 24 3/8 - 27 1/2"	62 - 70 cm 24 3/8 - 27 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	75 - 83 cm 29 1/2 - 32 5/8"	75 - 83 cm 29 1/2 - 32 5/8"



**SIZING:**

Measure circumference of the ankle, right above the round bones (malleoli) (cB), around the fullest part of the thigh (cG) and around the largest part of the hip (cH). Measure thigh height (IG).

SIZE	I	II	III	IV	V	VI	VII
cH	85 - 97 cm 33 1/2 - 38 1/8"	90 - 102 cm 35 3/8 - 40 1/8"	95 - 107 cm 37 3/8 - 42 1/8"	100 - 112 cm 39 3/8 - 44 1/8"	105 - 117 cm 41 3/8 - 46"	110 - 122 cm 43 1/4 - 48"	112 - 125 cm 44 1/8 - 49 1/4"
cG	45 - 53 cm 17 5/8 - 20 7/8"	48 - 56 cm 18 7/8 - 22"	51 - 59 cm 20 - 23 1/8"	54 - 63 cm 21 1/4 - 24 3/4"	57 - 66 cm 22 3/8 - 26"	60 - 70 cm 23 5/8 - 27 1/2"	63 - 73 cm 24 3/4 - 28 3/4"
cB	18 - 20 cm 7 - 7 7/8"	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"
IG	62 - 70 cm 24 3/8 - 27 1/2"	62 - 70 cm 24 3/8 - 27 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	75 - 83 cm 29 1/2 - 32 5/8"	75 - 83 cm 29 1/2 - 32 5/8"





2883

MATERNITY COMPRESSION PANTYHOSE

Indication

- Mild varicosity
- Onset of vein related problems
- Prevent edema and venous disorder during pregnancy




Size: I - VII

Color



Features:

- A. Mesh design in abdominal region won't stress the fetus.
- B. Button fastener on waist belt for easy adjustment.
- C. Reinforced stocking sole for long-wearing performance.



SIZING:

Measure circumference of the ankle, right above the round bones (malleoli) (cB) and around the fullest part of the thigh (cG). Measure thigh height (IG).

SIZE	I	II	III	IV	V	VI	VII
cG	45 - 53 cm 17 5/8 - 20 7/8"	48 - 56 cm 18 7/8 - 22"	51 - 59 cm 20 - 23 1/8"	54 - 63 cm 21 1/4 - 24 3/4"	57 - 66 cm 22 3/8 - 26"	60 - 70 cm 23 5/8 - 27 1/2"	63 - 73 cm 24 3/4 - 28 3/4"
cB	18 - 20 cm 7 - 7 7/8"	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"
IG	62 - 70 cm 24 3/8 - 27 1/2"	62 - 70 cm 24 3/8 - 27 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	75 - 83 cm 29 1/2 - 32 5/8"	75 - 83 cm 29 1/2 - 32 5/8"

Voyage Series

After long distance travel, some passengers may develop Economy Class Syndrome (Deep Vein Thrombosis-DVT). Wearing OPPO Voyage Compression Stockings and exercising your legs during the trip help improve blood circulation, relieve pain and edema of lower limbs, and counteract the sedentary conditions that can contribute to DVT. Sophisticated design for men and women respectively for better fit and a contoured design at the calf muscles ensure maximum comfort during a long Journey.





2826

TRAVEL STOCKINGS

Indication

- Improve blood circulation in lower limbs
- Prevent leg swelling, edema and discomfort
- Reduce the risk of DVT during long distance travel

Size: I - VII

Color (Female)

Color (Male)



2821

TRAVEL STOCKINGS

Indication

- Improve blood circulation in lower limbs
- Prevent leg swelling, edema and discomfort
- Reduce the risk of DVT during long distance travel

Size: I - VII

Color (Female)

Color (Male)



2800

TRAVEL STOCKINGS

Indication

- Improve blood circulation in lower limbs
- Prevent leg swelling, edema and discomfort
- Reduce the risk of DVT during long distance travel

Size: I - VII


Color (Female)

Color (Male)



SIZING:

Measure circumference of the ankle, right above the ankle bones (malleoli) (cB) and around the fullest part of the calf (cC). Measure knee height (ID).



SIZE	I	II	III	IV	V	VI	VII
Female							
cC	31 - 35 cm 12 1/4 - 13 3/4"	33 - 38 cm 13 - 15"	35 - 40 cm 13 3/4 - 15 3/4"	37 - 42 cm 14 1/2 - 16 1/2"	39 - 44 cm 15 1/4 - 17 1/4"	41 - 47 cm 16 1/8 - 18 1/2"	43 - 49 cm 16 7/8 - 19 1/4"
cB	18 - 20 cm 7 - 7 7/8"	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"
ID	35 - 39 cm 13 3/4 - 15 3/8"	35 - 39 cm 13 3/4 - 15 3/8"	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	40 - 45 cm 15 3/4 - 17 5/8"	40 - 45 cm 15 3/4 - 17 5/8"
Male							
cC	33 - 38 cm 13 - 15"	35 - 40 cm 13 3/4 - 15 3/4"	37 - 42 cm 14 1/2 - 16 1/2"	39 - 44 cm 15 1/4 - 17 1/4"	41 - 47 cm 16 1/8 - 18 1/2"	43 - 49 cm 16 7/8 - 19 1/4"	45 - 51 cm 17 5/8 - 20"
cB	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"	32 - 34 cm 12 1/2 - 13 3/8"
ID	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	42 - 47 cm 16 1/2 - 18 1/2"	42 - 47 cm 16 1/2 - 18 1/2"	42 - 47 cm 16 1/2 - 18 1/2"	43 - 48 cm 16 7/8 - 18 7/8"	43 - 48 cm 16 7/8 - 18 7/8"

## Eminence Series

OPPO Eminence Compression Stockings provide comfortable support for tired legs during work, exercise or walking, effectively reducing the risk of varicose veins and related symptoms. The classic ribbed design makes the stockings look like ordinary socks which are ideal for everyday wear both in business or casual shoes. The generous space for the toe and heel ensures non-restricting movement for superior comfort when worn over an extended period of time.







2828

KNEE HIGH COMPRESSION STOCKINGS

Indication

- Mild varicosity
- Leg pain, edema and muscle weakness

Size: I - VII

Color



2804

KNEE HIGH COMPRESSION STOCKINGS

Indication

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear

Size: I - VII

Color



2829

KNEE HIGH COMPRESSION STOCKINGS

Indication

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear

Size: I - VII

Color



2873

THIGH HIGH COMPRESSION STOCKINGS

Indication

- Mild varicosity
- Leg pain, edema and muscle weakness

Size: I - VII

Color



2863

THIGH HIGH COMPRESSION STOCKINGS

Indication

- Moderate to severe varicose veins
- Spider veins
- Post-surgical wear


Size: I - VII

Color



SIZING:


Measure circumference of the ankle, right above the ankle bones (malleoli) (cB) and around the fullest part of the calf (cC). Measure knee height (ID).



SIZE	I	II	III	IV	V	VI	VII
cC	33 - 38 cm 13 - 15"	35 - 40 cm 13 3/4 - 15 3/4"	37 - 42 cm 14 1/2 - 16 1/2"	39 - 44 cm 15 1/4 - 17 1/4"	41 - 47 cm 16 1/8 - 18 1/2"	43 - 49 cm 16 7/8 - 19 1/4"	45 - 51 cm 17 5/8 - 20"
cB	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"	32 - 34 cm 12 1/2 - 13 3/8"
ID	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	42 - 47 cm 16 1/2 - 18 1/2"	42 - 47 cm 16 1/2 - 18 1/2"	42 - 47 cm 16 1/2 - 18 1/2"	43 - 48 cm 16 7/8 - 18 7/8"	43 - 48 cm 16 7/8 - 18 7/8"

SIZING:

Measure circumference of the ankle, right above the ankle bones (malleoli) (cB), around the fullest part of the calf (cC) and around the fullest part of the thigh (cG). Measure thigh height (IG).



SIZE	I	II	III	IV	V	VI	VII
cG	51 - 59 cm 20 - 23 1/4"	54 - 62 cm 21 1/4 - 24 3/8"	57 - 66 cm 22 3/8 - 26"	60 - 69 cm 23 1/2 - 27 1/8"	63 - 73 cm 24 3/4 - 28 3/4"	66 - 76 cm 26 - 29 7/8"	69 - 79 cm 27 1/8 - 31"
cC	33 - 38 cm 13 - 15"	35 - 40 cm 13 3/4 - 15 3/4"	37 - 42 cm 14 1/2 - 16 1/2"	39 - 44 cm 15 1/4 - 17 1/4"	41 - 47 cm 16 1/8 - 18 1/2"	43 - 49 cm 16 7/8 - 19 1/4"	45 - 51 cm 17 5/8 - 20"
cB	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"	32 - 34 cm 12 1/2 - 13 3/8"
IG	70 - 80 cm 27 1/2 - 31 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	75 - 85 cm 29 1/2 - 33 1/2"	75 - 85 cm 29 1/2 - 33 1/2"	75 - 85 cm 29 1/2 - 33 1/2"	79 - 88 cm 31 - 34 1/2"	79 - 88 cm 31 - 34 1/2"





## Anti-embolism Series

OPPO Anti-embolism Compression Stockings provide equal pressure of 18 mmHg on the lower limbs. It is specially designed for bedridden patients during hospitalization or at home to help blood flow back to the heart and prevent the formation of blood clots in the vein, reducing the risk of thrombosis and resultant symptoms.



2805

KNEE HIGH ANTI-EMBOLISM STOCKINGS

Indication

- Prevent venous stasis
- Suitable for post-surgical or bedridden patients
- Reduce the formation of thrombosis or embolism

Size: I - VII

Color



2862

THIGH-HIGH ANTI-EMBOLISM STOCKINGS

Indication

- Prevent venous stasis
- Suitable for post-surgical or bedridden patients
- Reduce the formation of thrombosis or embolism


Size: I - VII

Color



Features:

- A.Comfortable wide band to prevent slipping and binding.
- B.Colored marker identifying the correct heel and foot area to provide constant compression.
- C.Opening in the ball-of-foot for easier monitoring of circulation without removing the stockings.




SIZING:

2805

Measure circumference of the ankle, right above the ankle bones (malleoli) (cB) and around the fullest part of the calf (cC). Measure knee height (ID).

2862

Measure circumference of the ankle, right above the ankle bones (malleoli) (cB), around the fullest part of the calf (cC) and around the fullest part of the thigh (cG). Measure thigh height (IG).



SIZE	I	II	III	IV	V	VI	VII
cG	45 - 53 cm 17 5/8 - 20 7/8"	48 - 56 cm 18 7/8 - 22"	51 - 59 cm 20 - 23 1/8"	54 - 63 cm 21 1/4 - 24 3/4"	57 - 66 cm 22 3/8 - 26"	60 - 70 cm 23 5/8 - 27 1/2"	63 - 73 cm 24 3/4 - 28 3/4"
cC	31 - 35 cm 12 1/4 - 13 3/4"	33 - 38 cm 13 - 15"	35 - 40 cm 13 3/4 - 15 3/4"	37 - 42 cm 14 1/2 - 16 1/2"	39 - 44 cm 15 1/4 - 17 1/4"	41 - 47 cm 16 1/8 - 18 1/2"	43 - 49 cm 16 7/8 - 19 1/4"
cB	18 - 20 cm 7 - 7 7/8"	20 - 22 cm 7 7/8 - 8 5/8"	22 - 24 cm 8 5/8 - 9 1/2"	24 - 26 cm 9 1/2 - 10 1/4"	26 - 28 cm 10 1/4 - 11"	28 - 30 cm 11 - 11 3/4"	30 - 32 cm 11 3/4 - 12 1/2"
IG	62 - 70 cm 24 3/8 - 27 1/2"	62 - 70 cm 24 3/8 - 27 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	70 - 80 cm 27 1/2 - 31 1/2"	75 - 83 cm 29 1/2 - 32 5/8"	75 - 83 cm 29 1/2 - 32 5/8"
ID	35 - 39 cm 13 3/4 - 15 3/8"	35 - 39 cm 13 3/4 - 15 3/8"	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	39 - 44 cm 15 3/8 - 17 1/4"	40 - 45 cm 15 3/4 - 17 5/8"	40 - 45 cm 15 3/4 - 17 5/8"



## APPENDIX

Specification		Caring Series																Voyage Series				Eminence Series				Anti-embolism Series	
		2830	2851	2871	2891	2883	2824	2852	2875	2880	2802	2840	2865	2882	2803	2861	2884	2826	2821	2800	2828	2873	2829	2804	2863	2805	2862
Compression level	I	●	●	●	●	●	○	○	○	○	○	○	○	○	○	○	○	●	○	○	○	●	○	○	○	18 mmHg	
	II	○	○	○	○	○	●	●	●	●	●	●	●	●	○	○	○	○	●	○	○	○	●	○	○		
	III	○	○	○	○	○	○	○	○	○	○	○	○	○	●	●	●	○	○	○	○	○	○	○	○		
Length	Knee high	●	○	○	○	○	●	○	○	○	○	○	○	○	○	○	○	●	○	○	○	○	○	○	○	○	○
	Mid thigh high	○	●	○	○	○	○	●	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Thigh high	○	○	●	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
Toe	Pantyhose	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Open	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Closed	●	●	●	●	●	●	●	●	●	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
Color	Beige	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●
	Black	○	●	●	○	○	●	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Gray	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
Indication	Prevent economy class syndrome	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Edema / Swollen legs	●	●	●	●	●	●	●	●	●	●	●	●	●	○	○	○	○	○	○	○	○	○	○	○	○	○
	Muscle weakness	●	●	●	●	●	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Leg pain	●	●	●	●	●	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Mild varicosity	●	●	●	●	●	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Pregnancy	●	●	●	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Post-surgical wear	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Moderate varicose veins	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Scleroderma	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Severe edema	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
Contraindication:	Recurrent varicose veins	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Severe varicosity during pregnancy	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
	Chronic venous insufficiency	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○	○
Contraindication:		Dermatitis, post-operative vein ligation, gangrene, recent skin graft, suspected pre-existing deep vein thrombosis, ischemic vascular disease, severe lower-extremity edema.																									



Support Your  
Healthy Lifestyle